

# ORDER INFORMATION



UK Distributors  
Hospital Lamp Supplies  
(A Division of Hybec Limited)  
Unit 9 Barrington Park  
Leycroft Road  
Leicester  
LE4 1ET

Tel. 0116 235 8818  
Fax. 0116 235 8810

e-mail hls@hybec.com

Website: www.hybec.com

**Price £325.00 ex VAT**

**UK carriage £9.99 ex VAT**

Prices correct at the time of publication



MD80938



# TECHNICAL DATA

## Technical Data

Dimensions (case) LxWxH: 290x390x110 mm  
Dimensions (Idromed 5) LxWxH: 200x170x87 mm  
Total weight 2.5kg (appliance, case and accessories)  
Power Supply: 230v / 50Hz  
UK 3 pin plug

Subject to technical alterations.

Supplied with a standard UK plug in 230v, English instructions, English language display.

We strongly recommend that you seek professional medical guidance before purchasing any dermatology products

NOTE: This product is VAT Exempt to most private users in the UK and Europe. However, in order to qualify for this you must complete the form supplied and declare that you are buying this equipment for your own (or spouse, child, etc) medical use.

As official UK distributors to hospitals and home users, we can offer the full 2 years manufacturer's warranty.

### In which cases is idromed® 5 applicable?

Besides the classic treatment of hyperhidrosis the idromed® 5 can also be used in the following cases:

- Gram-negative foot infections
- Dyshidrotic eczema
- Andrews syndrome (Pustulosis palmaris et plantaris)
- Acrocyanosis
- Arthrosis

### In which cases you may not use idromed® 5 ?

In the case of the following conditions, the idromed® 5 should not be used or should only be used after consulting with your physician:

- Arrhythmia
- Electronic implants (pacemaker)
- Metallic implants in the area of the current flow
- Metallic intrauterine implants (contraceptive coil)
- Pregnancy
- Extensive dermatological lesions

# HYPERHIDROSIS TREATMENT



**The drug-free iontophoresis therapy against excessive sweating**





**No risk of allergic reactions thanks to high-quality plate electrodes.**

The plate electrodes are made of aluminium. They are connected to the appliance by cables and are leading the direct current into the water-filled treatment basins.



**Practical handling**

Using the rotary control button and intuitive character display with built-in timer, a user can set up therapy without additional assistance



**All-in-one treatment basins and case**

The idromed® 5 is packed in a case, which not only helps to store and transport the appliance comfortably and safely - the two parts of the case can be separated and are used as treatment basins for the iontophoresis therapy.



**Standard Idromed 5 GS complete medication kit as shown above**

Economical and comfortable therapy which can equally be used in the home or in a Hospital or Medical Practice.

The amperage is increased until the patient feels a light 'tingling' sensation. The weak current slows down the sweat gland secretion without damaging the glands. Dependent upon individual circumstances, initially a daily 15 minute treatment is applied. Ususally 10 to 15 applications will normalize the moisture content of the skin, and thereafter it is sufficient to treat the affected area once or twice a week.

As with all dermatology products, it is recommended that guidance is sought from your doctor.

**Optional sponge electrodes for armpits**

With the special sponges it is also possible to treat hyperhidrosis in the armpits.



Telephone for prices of accessories

**Eligibility Declaration by a Disabled Person for products which are designated as eligible for VAT exemption as detailed and confirmed by Hybec Limited**

Please note there are penalties for making false declarations – and Hybec will only process VAT exemptions for people with permanent disabilities on and for products which qualify by law – and via manufacturer statement – as being solely designed for people with permanent disability. We will contact you immediately should we receive a declaration form not congruous to the requirements of the VAT exemption plan.

**CUSTOMER:**

If you are in any doubt as to whether you are eligible to receive goods or services zero-rated for VAT you should consult HM Revenue & Customs Notice 701/7 Reliefs for Disabled People or contact their National Advice Service on 0845 010 9000 before signing this declaration.

I (full name) \_\_\_\_\_

of (address) \_\_\_\_\_

**declare that**

• I am chronically sick or have a disabling condition by reason of (give full and specific description of your condition) :

• I am receiving from Hybec Limited (of Unit 9 Barrington Park, Leycroft Road, Leicester LE4 1ET) the following goods which are being supplied to me for my domestic or my personal use:

**Dr Hoenle Idromed 5**

I, the customer, further declare that I claim relief from value added tax (VAT)

Signature \_\_\_\_\_

Date: \_\_\_\_\_

**Customer, please return your signed and completed form by one of three ways:**

**Fax** +44 (0) 116 235 8810, **Post** (Hybec Limited, Unit 9 Barrington Park, Leycroft Road, Leicester LE4 1ET), or by **e-mailing a scanned pdf, jpg or gif formatted document** to enquiries@hybec.com.

**SUPPLIER**

I (full name): \_\_\_\_\_ of Hybec Limited, Unit 9 Barrington Park, Leycroft Road, Leicester LE4 1ET am supplying to the person named above, the following goods

**Dr Hoenle Idromed 5**

Signature: \_\_\_\_\_

Date \_\_\_\_\_